

PLASMA LIPIDS IN-VITRO FILTERING METHOD AND APPARSTUS

ABSTRACT OF THE DISCLOSURE

An in-vitro blood plasma lipids filtering method includes the following steps: separating out the blood plasma from the blood collection; flushing the apparatus carrying out the method with saline solution; controlling the temperature and pressure of the blood plasma; passing the blood plasma to screening procedure for filtering; and feeding the blood plasma back to the blood. The method is clearly effective and accurate, quick response indication, more secure and safer, more tolerant, and the treatment time is short.